




7300 Cessna Drive  
Greensboro, NC 27409  
Phone: 877.572.9985  
Fax: 877.572.9986  
www.clinicalchoice.com



## DISPOSABLE ENDOSCOPIC CLEANING BRUSHES

### Product Description

Clinical Choice™ Disposable Endoscopic Cleaning Brushes are available in 2 convenient styles to meet all of your cleaning needs.

Instrument channel brush tips are rounded to prevent endoscope channel damage and ensure smooth insertion. Brushes are available in large and small diameters with durable tubing providing excellent rigidity and maximum cleaning.

Packaged 50/box • Non-sterile • Single Use Only

### Features

- » Color coded for convenient, quick identification of catheter and/or brush size
- » Rounded brush tip prevents endoscope channel damage and provides maximum debris removal
- » Durable tubing and brush is easy to use and sturdy for maximum cleaning
- » Single use reduces risk of cross contamination.

### CBSE-2



Small diameter disposable cleaning brush (channel only)

Catheter Diameter	Catheter Length	Bristle Diameter
1.5mm	180cm	3mm

### Vendor Set Up Info

#### Federal ID Number

56-2245673

#### Website

www.clinicalchoice.com

#### Phone

877.572.9985

#### Fax

877.572.9986

#### Mailing Address

7300 Cessna Drive  
Suite A  
Greensboro, NC 27409

#### Remit To Address

PO Box 890368  
Charlotte, NC 28289-0368

### Place Your Order

#### Orders may be placed by:

##### Phone

877.572.9985

##### Fax

877.572.9986

##### E-Mail

orders@clinicalchoice.com



Contact Customer Service with any questions you may have.

### COMBO-1



Double-ended disposable cleaning brush (valve and channel combo)

Catheter Diameter	Catheter Length	Bristle Diameter
1.65mm	230cm	11mm & 5mm

CLINICAL CHOICE DISPOSABLE ENDOSCOPIC CLEANING BRUSHES

For Pricing Call 877-572-9985