

READY, SET, GO!

About Us

Clinical Choice[®] is a woman-owned business. Our number one goal is to provide you significant savings while maintaining a healthy environment for your patients. Our OEM quality compatible replacement filters and line of disposable products provide you with 30% - 50% savings over the OEM without compromising patient safety. We are pleased to have the opportunity to earn your business.

Implementation Made Easy

At Clinical Choice[®], we want to make it easy to enjoy immediate cost savings. With our **Ready, Set, Go** Program, you are provided a pre-assigned account number and access to all information you will need to quickly set us up in your system. We accept orders by phone, fax, or via email: orders@clinicalchoice.com or GHX. We offer same-day shipping and a 100% satisfaction guarantee. Our friendly and knowledgeable customer service staff is ready to answer any questions you have.



Compatible filters for OER PRO[®] System

SET

Vendor Set Up Info

Federal ID No.
56-2245673

Website
www.clinicalchoice.com

Phone
877.572.9985

Fax
877.572.9986

Mailing Address
5603 W. Friendly Ave.
Ste B. PMB 255
Greensboro, NC 27410

Remit to Address
PO Box 890368
Charlotte, NC 28289-0368

GO

Place Your Order

Orders may be placed by:

Phone
877.572.9985

Fax
877.572.9986

Email
orders@clinicalchoice.com



Contact Customer Service with any questions you may have.

FOR ADDITIONAL ASSISTANCE OR INFORMATION ON FILTER COMBINATION PACKAGES,
PLEASE CONTACT CUSTOMER SERVICE AT 877.572.9985

STERIS® Compatible

Type	OEM Model No.	CC Model No.	Description (All quantities are for single units unless otherwise noted)
Filter	MF01-0014	MFO-100	DSD and SSD 1 micron filter
Filter	CF01-0017	MFO-150	1 micron pre-filter for (5000, MV, CLM, CER)
Filter	MF01-0015	MFO-200	DSD and SSD 0.45 micron filter
Filter	MF01-0003	MFO-250	DSD 5 micron inline pre-sediment water filter
Filter	MF01-0011	MFO-300	DSD and SSD disinfectant filter (12/pkg)
Filter	MF01-0028	MFO-375	DSD and SSD air filter (8/pkg)
Filter	MF01-0013	MFO-500	DSD and SSD 0.2 micron bacterial retention grade filter
Filter	CF01-0018	MFO-600	0.2 micron post filter for (5000, MV, CLM, CER)
Filter	MF01-0078	MFO-750	DSD EDGE complete filter package: (1) MFO-375, (4) MFO-100 (4) MFO-200, (2) MFO-500
Filter	MF01-0070	ADP-100	1 micron Nominal pre-filter
Filter	MF01-0071	ADP-200	0.4 micron Absolute pre-filter
Filter	MF01-0069	ADP-500	0.1 micron Absolute pre-filter

OLYMPUS® OER PRO® Compatible

Type	OEM Model No.	CC Model No.	Description (All quantities are for single units unless otherwise noted)
Filter	MF01-0014PL	OLY-100	1 micron pre-sediment filter
Filter	MF01-0015PL	OLY-200	0.45 micron pre-filter
Filter	MAJ-823	OLY-300	Internal air filter
Filter	MAJ-824	OLY-400	Dual layer bacterial retentive filter with centering ring

STERIS® Compatible

Type	OEM Model No.	CC Model No.	Description (All quantities are for single units unless otherwise noted)
Filter	A1501	MF-STE A	"A" pre-filter for Reliance Processing System
Filter	A1504	MF-STE B	"B" pre-filter
Filter	A1551	MF-STE A/B	"A" and "B" combination pre-filter set
Filter	A1501E	MF-STE A	"A" pre-filter for STERIS System 1E™
Filter	A1562	MF-STE B1E	"B" pre-filter for STERIS System 1E™
Filter	A1567	MF-STE A/B1E	"A" and "B" combination pre-filter set for STERIS System 1E™
Filter	A1503	MF-STE D	Internal air filter

CUSTOM ULTRASONICS® Compatible

Type	OEM Model No.	CC Model No.	Description (All quantities are for single units unless otherwise noted)
Filter	16030	CU-50	Water and/ or air adapter insert (50/pkg)
Filter	13015	CU-100	25 micron pre-sediment filter (6/pkg)
Filter	13016	CU-200	5 micron pre-sediment filter (6/pkg)
Filter	13017	CU-300	0.1 micron bacteria filter
Filter	21017	CU-400	20" cartridge carbon block filter

ASP® AER Compatible

Type	OEM Model No.	CC Model No.	Description (All quantities are for single units unless otherwise noted)
Filter	20303	MFJ-100	1 micron pre-sediment filter (30/pkg)
Filter	20304	MFJ-200	0.2 micron bacteria filter

Disclaimer: The presentation of such product names, product descriptions, and marks herein does not represent any endorsement or approval of Clinical Choice® or its products or services by the owners of such product names, product descriptions, or marks, or any agreement or affiliation between Clinical Choice® and such owners and are intended to be informational in nature.