

# BiteBlock

Top-rated solutions that safeguard your patients and endoscopy equipment.



## Adult BiteBlock

ORDER NUMBER

**CCBB-L**

SIZE

**20x27 mm - 60 Fr**

QUANTITY

**Sold in 50/bag**



## Pediatric BiteBlock

ORDER NUMBER

**CCBB-S**

SIZE

**14x20 mm - 42 Fr**

QUANTITY

**Sold in 50/bag**

READY TO PLACE AN ORDER?

E: [orders@clinicalchoice.com](mailto:orders@clinicalchoice.com)

T: (877) 572-9985


F: (877) 572-9986

GOT QUESTIONS?

You may reach our Customer Service via phone or [info@clinicalchoice.com](mailto:info@clinicalchoice.com)



Visit [clinicalchoice.com](http://clinicalchoice.com) for information on all product offerings: GI, SPD/CS and OR infection control products.



## Why choose Clinical Choice BiteBlocks?



### Seamless Form and Function

Our Bite Blocks reduce endoscope friction, provide convenient access for suctioning, and allow easy adjustments throughout procedures.



### Shaped for Comfort

Our Bite Blocks are flexible and adaptable. They conform to the sides of the mouth and face, allowing them to maintain a proper position during endoscopy procedures.



### Designed for Safety

Our Bite Blocks and adjustable straps are latex free. The dental retention rim reduces the risk of gum damage and prevents patients from biting down on the endoscope insertion tube.

## The Clinical Choice Advantage

For over 21 years, we have offered significant savings, delivered high-quality products and provided excellent customer support. We currently serve over 1,200 hospitals, endoscopy and ambulatory surgery centers. We are truly proud to state that our products are being used by thousands of professionals daily.